

High-quality consumables



The performance of an analytical system is only as good as its consumables. That is why our consumables have been specifically developed to complement our analytical systems. Because precise analytical results can only be guaranteed if consumables and systems are designed to work together.

Kjeldahl analysis

Consumables for your daily work:



KJELCAT tablets

- + Immediate availability due to production in Germany
- + Stably pressed tablets for immediate use - without weighing
- + Long life time
- + Convenient packaging for easy removal
- + No pulverisation, no waste
- + Visual inspection of the alkalinisation of the sample through colour change
- + Wide range of KJELCAT variants, tailored to your application

Order no.	Type	Composition
12-0326	Se	3,5 g K_2SO_4 + 0,0035 g Se
12-0328	Cu	5,0 g K_2SO_4 + 0,5 g $CuSO_4 \times 5H_2O$
12-0329	CuTi	5,0 g K_2SO_4 + 0,15 g $CuSO_4 \times 5H_2O$ + 0,15 g TiO_2
12-0336	Cu light	5,0 g K_2SO_4 + 0,1 g $CuSO_4 \times 5H_2O$
12-0334	CK	3,5g K_2SO_4 + 0,4g $CuSO_4 \times 5H_2O$



Antifoam tablets

- + Effective prevention of foaming
- + Easy dosing thanks to tablet form

Order number 12-0332



Weighing paper


- + Nitrogen-free paper
- + Easy to remove and shape: makes sample weighing and transfer easier
- + Residue-free decomposition during oxidation

Order number 12-0353



Wide range of glass ware

- ✦ Particularly good sealing due to robust ground joint
- ✦ Robust quartz glass with extra thick glass walls for a long life time of digestion tubes
- ✦ Versatile glass variations tailored to the demands of your application



Buffer solution

- ✦ Standardised and ready-to-use pH buffer solutions for calibrating the pH electrode

Order number 1000176 pH-4
Order number 1000177 pH-7



Electrolyte solution KCl

- ✦ Ready-to-use for an optimal storage life of pH electrode

Order number 1000178



Activated carbon for refilling the water-free gas scrubber VACUSOG

- ✦ The kit ensures the exact filling quantity and a tight sealing for the activated carbon trap
- ✦ The evenly applied KOH coating achieves neutralisation of up to 1,000 samples
- ✦ Reliable neutralisation of acidic vapours for more safety in your laboratory

Order number 12-0785



Consumables by CG

- Top-quality analyses
- High-quality reagents
- Guaranteed performance with C. Gerhardt systems
- In accordance with ISO 9001
- Proven development process
- In-house quality assurance



Reliable partner

- Innovative products
- Many years of experience
- Personalised consultation
- Training for your employees
- Attractive quantity discounts



Our promise to you

- Smooth order process
- Secure supply chain
- Reliable availability
- Flexible service concepts
- Made in Europe
- Excellent price-performance ratio

You can also find all of our consumables on our website. Conveniently put together a selection of consumables for your daily work here.



www.gerhardt.de/en/products/accessories-consumables/

Contact:
lead@gerhardt.de

or your authorized local distributor

www.gerhardt.de/en