

The Kjeldahl method:

REFERENCE METHOD FOR NITROGEN AND PROTEIN ANALYSIS

Analytical process for the determination of nitrogen and protein according to Jean Kjeldahl in national and international standards and methods



Meat

- ✦ ISO 937
- ✦ AOAC 928.08
- ✦ AOAC 981.10
- ✦ LFGB § 64 method L 06.00-7



Beer

- ✦ AOAC 920.53
- ✦ AOAC 945.23
- ✦ AOAC 950.09
- ✦ AOAC 950.10
- ✦ MEBAK 1.5.2.1
- ✦ MEBAK 3.1.4.5.1.1
- ✦ MEBAK 3.1.4.5.2.1



Agriculture / Fertilizer

- ✦ AOAC 955.04
- ✦ AOAC 959.04
- ✦ AOAC 970.02
- ✦ AOAC 978.04



Environmental

- ✦ ISO 5663
- ✦ ISO 5664
- ✦ DIN 38406-5
- ✦ AOAC 973.48
- ✦ EPA 351.3
- ✦ ASTM D1426
- ✦ ASTM D1783



Milk and dairy products

- ✦ ISO 8968-1
- ✦ ISO 8968-3
- ✦ ISO 8968-4
- ✦ AOAC 930.29
- ✦ AOAC 939.02
- ✦ AOAC 991.20 - 23
- ✦ LFGB § 64 method L 01.00-10/1



Food

- ✦ ISO 1871
- ✦ ISO 3188
- ✦ ISO 20483
- ✦ AOAC 920.152
- ✦ AOAC 925.31
- ✦ AOAC 930.33
- ✦ AOAC 932.08
- ✦ AOAC 940.25
- ✦ AOAC 950.48



Grain

- ✦ AOAC 979.09
- ✦ ICC Standard No. 105/2



Rubber

- ✦ ISO 1656
- ✦ ASTM D3533



Industry

- ✦ ISO 3332
- ✦ ASTM D3228
- ✦ ASTM UOP384-76



Feedstuff

- ✦ ISO 5983-2
- ✦ EC No 152/2009
- ✦ AOAC 954.01
- ✦ AOAC 976.05
- ✦ AOAC 984.13
- ✦ AOAC 988.05
- ✦ AOAC 2001.11
- ✦ GAFTA No. 130



Your Kjeldahl method is not listed here? Please contact us: C. Gerhardt strives to keep these standards and methods up to date. Together we optimize your analytical procedure and automate it with our systems.